Chapter 23.2 - Ester Lab

<u>Part I</u>

Procedure A

- Obtain and label a micro-tube A. Make sure the tube is clean and dry.
- 2. Using a clean pipet, *place 10 drops* of *isopentyl alcohol* into the micro-tube.
- 3. Using a clean pipet, place 10 drops of glacial acetic acid into the micro-tube.
- 4. Add *one drop* of concentrated *sulfuric acid* into the micro-tube.
- 5. Close the micro-tube and set aside.

Procedure B

- 1. Obtain and label a micro-tube B. Make sure the tube is clean and dry.
- 2. Using a clean pipet, place 10 drops of octanol alcohol into the micro-tube.
- 3. Using a clean pipet, place 10 drops of glacial acetic acid into the micro-tube.
- 4. Add *one drop* of concentrated *sulfuric acid* into the micro-tube.
- 5. Close the micro-tube and set aside.

Procedure C

- 1. Obtain and label a micro-tube C. Make sure the tube is clean and dry.
- 2. Using a scale, measure **0.15** grams of salicylic acid and add it to the micro-tube.
- 3. Using a clean pipet, place 12 drops of methyl alcohol to the micro-tube.
- 4. Add three drops of concentrated sulfuric acid into the micro-tube.
- Close the micro-tube and set aside.

Procedure D

- Obtain and label a micro-tube D. Make sure the tube is clean and dry.
- Using a clean pipet, place 10 drops of propanol alcohol (n-Propyl) into the micro-tube.
- 3. Using a clean pipet, place 10 drops of glacial acetic acid to the micro-tube.
- 4. Add *one drop* of concentrated *sulfuric acid* into the micro-tube.
- 5. Close the micro-tube and set aside.

Part II

- 1. Heat 50 ml of water in a beaker, using a hot plate, to 80 degrees Celsius.
- 2. Using tongs pour the hot water into a Styrofoam cup.
- 3. Place all four micro-tubes (#A D) with the top lid closed into cup of hot water.
- 4. After ten minutes, open micro-tubes A, B and D ONLY and carefully waft the odor.
- 5. For micro-tube <u>C</u>, open and add <u>3 drops of water</u> to the tube and carefully waft the odor.

23.2 – Ester Lab Report

| Micro-tube Letter | Name of Alcohol | Name of Acid | Odor of Ester |
|-------------------|-----------------|--------------|---------------|
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

| Questions | | | |
|---------------------------------------|----------------------------|---|-----------------------|
| acid being changed t | o -ate. For example, este | that is used to produce it w r ethyl benzoate is made fi the four esters that you pro | rom the alcohol ethyi |
| a.) | | | |
| b.) | | | |
| c.) | | | |
| d.) | | | |
| 2.) What substituted | hydrocarbon is found in a | all alcohois? | |
| 3.) Organic acids, like substitution? | the glacial acetic acid us | ed in this lab contain what | hydrocarbon |
| 4.) All esters contain | what hydrocarbon substit | tuted group? | |
| 5.) Give two example a.) | es of how esters are used | in our everyday lives. | |
| | | | |

b.)